

# Extraction of Tamoxifen and Metabolites from Urine and Plasma using Supported Liquid Extraction (SLE) prior to LC-MS/MS Analysis

Lee Williams, Rhys Jones, Helen Lodder, Geoff Davies, Steve Jordan, Richard Calverley, Claire Desbrow & Gary Douthwaite  
Biotage GB Limited, Distribution Way, Dyffryn Business Park, Ystrad Mynach, Hengoed, CF82 7TS, UK.

## Introduction

Tamoxifen is an important estrogen receptor antagonist primarily used in breast cancer therapy. More recently its action has also been shown to inhibit prostate cancer. Its mode of action also reduces the secondary effects linked to adsorption of androgen anabolic steroids. As a result the International Olympic Committee designated tamoxifen as a 'banned substance'. This widespread use along with various forms of misuse has led to the necessity of rapid and reliable methods for its analysis and quantification. Here we demonstrate a rapid and reliable 96-well Supported Liquid Extraction assay for the extraction of tamoxifen and metabolites from various human biological fluids.

## Experimental Procedure

### Reagents


Formic acid, ammonium hydroxide and tamoxifen were purchased from Sigma Chemical Co. (Poole, UK). 4-hydroxytamoxifen, N-desmethyltamoxifen and 4-hydroxy-N-desmethyltamoxifen (endoxifen) were purchased from Toronto Research Chemicals (North York, ON, Canada). Human plasma was obtained through the Welsh Blood Service (Pontyclun, UK). Urine was obtained from a healthy human volunteer. All solvents were HPLC grade from Fisher Scientific (Loughborough, UK).

### Sample Preparation

#### Supported Liquid Extraction Procedure

**Plate:** ISOLUTE SLE+ 200 Supported Liquid Extraction Plate, part number 820-0200-P01

**Sample pre-treatment:** Plasma or urine (100 µL) was diluted (1:1 (v/v)) with either 1% formic acid aq, 0.1% formic acid aq, H<sub>2</sub>O or 0.5M NH<sub>4</sub>OH. Spike concentrations were 200 ng/mL of plasma.



**Sample Application:** The pre-treated plasma/urine (200 µL) was loaded onto the plate, a pulse of vacuum applied to initiate flow and the samples left to absorb for 5 minutes.

**Analyte Elution:** Addition of 1 mL of various water immiscible extraction solvents. The extraction solvents evaluated were Hexane, MTBE, EtOAc, DCM.

**Post Extraction:** The eluate was evaporated to dryness and reconstituted in 500 µL of 0.1% formic acid 50:50 (v/v) H<sub>2</sub>O/MeOH prior to analysis.

### HPLC Conditions

**Instrument:** Waters Acquity UPLC (Waters Assoc., Milford, MA, USA).

**Column:** Acquity UPLC BEH C18 column (1.7µ, 100 x 2.1 mm id) (Waters Assoc., Milford, MA, USA).

**Mobile Phase:** 0.1% formic acid aq and 0.1% formic acid/MeOH at a flow rate of 0.5 mL/min.

**Gradient:** Isocratic 10%, 0.1% (v/v) formic acid aq and 90% 0.1% formic acid/MeOH.

**Injection Volume:** 5 µL.

**Column Temperature:** 35 °C.

### Mass Spectrometry

**Instrument:** Quattro Premier XE triple quadrupole mass spectrometer (Waters Assoc., Manchester, UK) equipped with an electrospray interface for mass analysis. Positive ions were acquired in the multiple reaction monitoring (MRM) mode.

**Desolvation Temperature:** 450 °C

**Ion Source Temperature:** 150 °C

**Collision Gas Pressure:** 3.56 x 10<sup>-3</sup> mbar

Scan Function	Analyte	Transition	Cone Voltage V	Collision Energy eV
1	Endoxifen	374.2 > 58.0	38	22
	4-OH-Tamoxifen	388.2 > 72.0	42	23
2	N-desmethyltamoxifen	358.2 > 58.0	37	21
	Tamoxifen	372.2 > 72.0	36	25

Table 1. Quattro Premier XE mass spectrometer parameters



## Results

Figures 1-4. demonstrate recoveries for the extraction of tamoxifen and metabolites from plasma. Recoveries > 70% were observed using a number of pre-treatment/extraction solvent protocols.




Figure 1. Tamoxifen recovery profile from plasma using 1% formic acid pre-treatment.




Figure 2. Tamoxifen recovery profile from plasma using 0.1% formic acid pre-treatment.




Figure 3. Tamoxifen recovery profile from plasma using H2O pre-treatment.




Figure 4. Tamoxifen recovery profile from plasma using 0.5M NH<sub>4</sub>OH pre-treatment.

Figures 5-8. demonstrate recoveries for the extraction of tamoxifen and metabolites from urine. Once again a variety of extraction protocols resulted in recoveries > 70% with corresponding RSDs < 10%.




Figure 5. Tamoxifen recovery profile from urine using 1% formic acid pre-treatment




Figure 6. Tamoxifen recovery profile from urine using 0.1% formic acid pre-treatment.






Figure 7. Tamoxifen recovery profile from urine using H<sub>2</sub>O pre-treatment.




Figure 8. Tamoxifen recovery profile from urine using 0.5M NH<sub>4</sub>OH pre-treatment.

## Conclusions

1. High reproducible recoveries were observed using various extraction combinations.
2. Hexane showed high recoveries of the less polar tamoxifen and desmethyltamoxifen but very low extraction efficiencies for the more polar metabolites.
3. Higher extraction efficiencies were observed for tamoxifen and desmethyltamoxifen from urine, suggesting that protein binding may occur in plasma.
4. EtOAc alone demonstrated acceptable extraction efficiencies from urine using 1% formic acid pH pre-treatment however, at higher loading pHs multiple solvents delivered recoveries > 80%.
5. Plasma extracts exhibited higher efficiencies using either H<sub>2</sub>O or 0.5M NH<sub>4</sub>OH pre-treatment.
6. This poster shows the applicability of SLE+ to extract tamoxifen and metabolites from plasma and urine.
7. RSDs were below 10%, demonstrating reproducible extraction and recoveries.

[www.biotage.com](http://www.biotage.com)

### United States and Canada

Tel: +1 704 654 4900  
Toll-Free: +1 800 446 4752  
ordermailbox@biotage.com

### United Kingdom, EIRE

Biotage  
Tel: +44 1992 501535  
order@eu.biotage.com

### Sweden

Biotage  
Tel: +46 18 56 59 00  
order@eu.biotage.com

### Japan

Biotage  
Tel: +81 422 281233  
order@biotage.co.jp

Copyright © 2010 Biotage. All rights reserved. All brand and product names are trademarks or registered trademarks of their respective companies. The information contained in this document is subject to change without notice.

