

Comparison of Biotage and Competitor Oligonucleotide Synthesis Methodologies



Synthesizing oligonucleotides often leads to the production of impurities that can be detrimental to downstream applications. These impurities are typically shorter, truncated versions of the full-length oligonucleotide. To ensure high-quality results, the solid phase synthesis and chromatographic purification processes must be carefully optimized and precisely controlled to minimize the presence of truncated side-products, leaving the oligonucleotide in a highly pure form. Biotage recognizes the need for continuous optimization of our synthesis and purification processes to ensure that our oligonucleotides are of the highest quality. This application note compares a Biotage oligonucleotide to the same sequence supplied by another vendor.

The Problem with Impurities

When synthesizing an oligonucleotide each iteration of the synthesis cycle can result in a small amount of incomplete coupling which leads to the formation of shorter failure sequences as by-products, i.e., oligonucleotide missing one or more nucleotides, so N -1 nucleotide, -2 nucleotides, -3 nucleotides and so on. If these failure sequences are not removed their presence can have negative impacts on downstream applications such as PCR or enzymatic reactions,

interfering with the target molecule's performance and distorting results. Furthermore, these impurities possess intrinsically similar chemical and physical properties to the desired full-length sequence and hence achieving chromatographic separation on a preparative scale is highly challenging. As such, in order to obtain high purity oligonucleotides a two-fold strategy is required: firstly, it is essential that all parameters associated with oligonucleotide synthesis are carefully controlled to ensure maximum coupling efficiency and minimal formation of side-products; secondly the purification process must be carefully optimized for each individual oligonucleotide sequence in order to ensure maximum resolution of any impurities which do form.

High Quality Oligonucleotide Preparation Techniques

At Biotage, we use our specialist chemical knowledge to ensure maximum efficiency during each step of the oligonucleotide synthesis cycle, with minimal side-product formation. In addition, we have developed a range of proprietary purification techniques designed specifically for removing failure sequences and other common impurities from synthesized oligos and delivering the highest levels of purity possible. Our expertise in purification ensures that our services consistently deliver extremely high-quality molecules at an affordable price point for our customers. Furthermore, investigations into materials from competitor companies have shown that our techniques are far superior when it comes to producing pure molecules free from unwanted side-products or contaminants.

Figures 1 and 2 show comparative capillary electrophoresis (CE) traces of the same oligonucleotide sequence – a 74-mer oligonucleotide containing standard phosphorothioate and 5'-phosphate modifications – synthesized by Biotage and also by a competitor: the presence of a series of deletion sequences is evident in the electropherogram of the competition sample. Final purities of the two samples as delivered to the customer were >95% for the Biotage material and >65% for the competitor.

Conclusion

In conclusion, utilizing proper purification techniques is essential for ensuring high-quality synthesized oligonucleotides free from any unwanted by-products or contaminants.

At Biotage, we leverage our expertise in purification technologies to provide highly pure molecules that you can confidently apply to your research and development projects without worry about potential problems caused by unwanted impurities in your sample. Investigate our range of products today and see how we can help you achieve success with your next project.




Figure 1. CE trace of an oligonucleotide synthesized by Biotage at >95% purity. The gel used for the CE analysis was ssDNA 100r kit and the length of the capillary was 30.2 cm. The voltage of separation that was be 9.0 kV over 60 min.




Figure 2. CE trace of the same oligonucleotide synthesized by a competitor at >65% purity. The gel used for the CE analysis was ssDNA 100r kit and the length of the capillary was 30.2 cm. The voltage of separation that was be 9.0 kV over 60 min.

EUROPE

Main Office: +46 18 565900
 Fax: +46 18 591922
 Order Tel: +46 18 565710
 Order Fax: +46 18 565705
 order@biotage.com
 Support Tel: +46 18 56 59 11
 Support Fax: + 46 18 56 57 11
 eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
 Toll Free: +1 800 446 4752
 Fax: +1 704 654 4917
 Order Tel: +1 800 446 4752
 Order Fax: +1 704 654 4917
 ordermailbox@biotage.com
 Support Tel: +1 800 446 4752
 us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
 Fax: +81 3 5627 3121
 jp_order@biotage.com
 jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 68162810
 Fax: +86 21 68162829
 cn_order@biotage.com
 cn-1-pointsupport@biotage.com

KOREA

Tel: +82 31 706 8500
 Fax: +82 31 706 8510
 korea_info@biotage.com
 kr-1-pointsupport@biotage.com

INDIA

Tel: +91 11 45653772
 india@biotage.com

Distributors in other regions
 are listed on www.biotage.com

Literature Number: PPS701

© 2023 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.