Biotage® V-10 Touch

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50 kg/110 lbs、CAROUSEL付き 55kg

Part Numbers

V10-2XX V-10 Touch System
V10-2SC V-10 Touch System with Solvent Manager and Automation

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
| Volume & Solvent | Method | Evaporation Time |
|-----------------------------------|----------------|------------------|
| 20 mL Vials | | |
| 8 mL NMP | Very-High Boil | 18 minutes |
| 8 mL DMSO | Very-High Boil | 15 minutes |
| 8 mL DMF | Hi-Boil | 4 minutes |
| 8 mL Water | Aqueous | 8 minutes |
| 8 mL Methanol | Volatile | 3.5 minutes |
| 8 mL DCM | High-Volatile | 4 minutes |
| 8 mL Ethanol | Volatile | 5 minutes |
| 8 mL 25% Piperidine in DMF | mxd Vol + HBP | 7 minute |

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| Part No. | 説明 |
|---------------------|---|
| エバポレーションシステム | |
| V10-2XX | V-10 Touch System |
| V10-2SC | V-10 Touch System with Solvent Manager and Automation |
| 設定のアップグレード | |
| 414865 | V-10 Solvent Manager Upgrade |
| 410876 | V-10 Automation Upgrade |
| 411816 | Gilson Liquid Handler GX271 - V-10 |
| アクセサリ | |
| 416009SP | Hi-vac, oil-free scroll pump, 100-240V |
| 413146 | Vacuum pump External Connection Kit |
| 411181 | Carousel Vial Holder 30/20 mL Vials |
| 411182 | Carousel Vial Holder 4 mL A-type HPLC Vials 14.5 mm Diameter (TYPE A/B) |
| 411183 | Carousel Vial Holder 4 mL Vials 16.5 mm Diameter (TYPE C) |
| 411985 | Carousel Vial Holder 16 mL Tall Vials 20.7 mm Diameter |
| 410647 | 4 mL Vial Upper Adapter Kit 5 pcs (Type A), 8.5 mm Neck Inner Diameter |
| 411031 | 4 mL Vial Lower Adapter Kit 2 pcs (Type A), 14.8 mm Vial Diameter |
| 410648 | 4mL Vial Upper Adapter Kit 5 pcs (Type B), 11.1 mm Neck Inner Diameter |
| 411032 | 4mL Vial Lower Adapter Kit 2 pcs (Type B), 14.1 mm Vial Diameter |
| 411033 | 4mL Vial Lower Adapter Kit 2 pcs (Type C), 16.3 mm Vial Diameter |
| 410649 | 4mL Vial Upper Adapter Kit 5 pcs (Type C), 12.0 mm Neck Inner Diameter |
| 412867 | 16 mL Upper Head Adapter Assembly |
| 412000 | 16 mL Vial Lower Adapter 2 pcs, (Type E), 20,7 mm Diameter |
| 410790 | Kit, Upper Vial Face Seal, 10 Pack |
| 412866 | Kit, Face Seal 16 mL, 10 pcs |

Service

| | |
|--------------|---|
| SER-V10T-SAP | Service Agreement - Priority - V-10 Touch Evaporator |
| SER-V10T-IN | Installation and Familiarization V-10 Touch Evap system |

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E-mail: Japan_info@biotage.com

Part Number: PPS425JPN_231129OSK

© 2020 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.

© Biotage 2020