


A fast heterocyclic three component synthesis of imidazo[1,2-a] annulated pyridines, pyrazines, pyrimidines and thiazoles under microwave conditions


Jacob Westman, Anders Franzén
Personal Chemistry, Kungsgatan 76, 753 18 Uppsala, Sweden

Introduction: Multicomponent reactions are becoming more popular nowadays since they provide a possibility to introduce a large set of diversity to different substance classes in one step. The imidazo[1,2-a]pyridazine, -pyridine, -pyrimidine and -thiazole structural moieties can be found in pharmacologically active compounds, such as anti-inflammatory agents, (Miroprofen)¹, calcium channel blockers² and antibacterials³. In 1998, three different research groups discovered a novel reaction at almost the same time⁴. They found (partly by serendipity) that when using 2-aminopyridines as the amine component together with isonitriles and aldehydes in the Ugi 4 component reaction (4CR) the formation of imidazo-pyridines was observed. The acid, the fourth component of the Ugi 4CR, is required as a catalyst. Under conventional conditions, the reaction takes 18-72 h depending on the method used to give the products in high yield.




In 1999 Varma *et al.* published this reaction under microwave and solvent-free conditions in a two-step fashion⁵. The two-step procedure introduces a limitation since alkyl iminium ions can be unstable and the solvent free procedure is inconvenient when it comes to library production. Especially when automated liquid handlers are used, but also due to the fact that many products are highly crystalline, making the isolation problematic.

We therefore followed the papers developed by Groebke *et al.* and Bienaymé *et al.* in order to find a solution phase protocol under microwave conditions and with short reaction times⁶. First we tried the method by Groebke *et al.* who used 2 eq. of acetic acid as catalyst. After up to 30 min at 170 °C in EtOH we only found traces of product (LC/MS analysis). Bienaymé used perchloric acid as the protic acid needed for the reaction. Perchloric acid may cause an explosion when heated and is therefore considered not to be very suitable for microwave heating. But in a closed pressurized vial and in only catalytic amounts we did not find any problems. 170 °C for 5 min was found to be a fairly general protocol.



The table below shows some structural variations tolerated by these reaction conditions. Both aromatic, aliphatic (even sterically hindered) and heteroaromatic aldehydes were used and gave the corresponding products in acceptable to good yields. The same is true for the isonitriles and heteroaromatic amidines used. It should be noted, though, that the yields presented in the table are based on the product peak area as compared to total peak area of the LC chromatogram (LC purity). As described above some of the products are highly crystalline, and in those cases were the products precipitated directly from the reaction mixture, the products were isolated by a simple filtration. No further purification was needed. (No recovery from the mother liquor implies that the total yield could be higher). When scaling up the synthesis [18 times], substrates 1b, 2j & 3e gave the isolated product in 70 % yield, as in the original protocol, see table below.



| Compound | Yield ^a | Compound | Yield ^a | Compound | Yield ^a | Compound | Yield ^a | Compound | Yield ^a |
|----------|--------------------|----------|--------------------|----------|--------------------|----------|--------------------|----------|--------------------|
| 1a2b3a | 83 | 1a2k3c | 97 | 1b2g3a | 50 | 1b2m3a | 83 | 1e2a3d | 63 (42) |
| 1a2f3a | 60 | 1a2l3a | 82 | 1b2l3a | 76 (71) | 1b2j3c | -- (70) | 1e2a3f | 62 (17) |
| 1a2f3b | 57 | 1a2l3c | 97 | 1b2l3c | 23 (58) | 1b2n3a | 95 | 1e2d3d | 64 (44) |
| 1a2g3a | 69 | 1a2m3a | 78 | 1b2l3e | 38 | 1b2n3c | 98 | 1e2d3f | 57 (9) |
| 1a2g3c | 69 (37) | 1a2m3b | 74 | 1b2j3a | 28 (54) | 1c2a3d | 39 | 1e2e3f | 49 |
| 1a2l3a | 75 | 1a2n3c | 95 | 1b2j3e | -- (71) | 1c2d3d | 42 (21) | 1e2h3d | 39 |
| 1a2l3c | 93 | 1b2b3a | 75 | 1b2k3a | 57 (51) | 1c2d3f | 35 | 1e2h3f | 81 (9) |
| 1a2j3a | 82 | 1b2c3e | 76 | 1b2k3c | 16 (66) | 1d2c3a | 62 | 1e2o3d | 46 |
| 1a2j3b | 29 (62) | 1b2c3a | 62 | 1b2k3e | 11 (76) | 1d2d3d | 42 | 1e2o3f | 55 |
| 1a2j3c | 27 (66) | 1b2g3e | 59 (53) | 1b2l3a | 89 (32) | 1d2e3d | 36 | | |
| 1a2k3a | 83 | 1b2f3a | 75 (46) | 1b2l3c | 56 (58) | 1d2h3d | 20 (9) | | |
| 1a2k3b | 76 (54) | 1b2f3c | 8 (60) | 1b2l3e | 56 | 1d2o3d | 39 | | |

^a = LC purity, in parentheses: isolated yield after filtration of the precipitated product.

References

- Maruyama, Y.; Anami, K.; Terasawa, M.; Goto, K.; Imayoshi, T.; Kadobe, Y.; Mizushima, Y. *Arzneimittel-Forsch.* 1981, 31, 1111-1118.
- Sanfilippo, P. J.; Urbanski, M.; Press, J. B.; Dubinsky, B.; Moore, J. B. *J. Med. Chem.* 1988, 31, 2221-2227.
- Rival, Y.; Grassy, G.; Michel, G. *Chem.Pharm.Bull.* 1992, 40, 1170-1176.
- Bienaymé, H.; Bouzid, K. *Angew. Chem. Int. Ed. Engl.* 1998, 37, 2234-2237.
- Varma, R. S.; Kumar, D. *Tetrahedron Lett.* 1999, 40, 7665-7669.
- Groebke, K.; Weber, L.; Mehlin, F. *Synlett* 1998, 661-663.; Bienaymé, H.; Bouzid, K. *Angew. Chem. Int. Ed. Engl.* 1998, 37, 2234-2237.