

# Reduction of Matrix Effects due to Phospholipids: A Comparison of Sample Preparation Techniques

**Lee Williams, Scott Merriman, Helen Lodder, Anne Howells, Steve Jordan, Claire Desbrow, Matthew Cleeve, Richard Calverley**

Argonaut Technologies, now a Biotage company, Dyffryn Industrial Estate, Ystrad Mynach, Mid Glamorgan, CF82 7RJ, UK

## Introduction

When analyzing drugs in biological fluids using LC-MS-MS, matrix effects due to endogenous phospholipids are commonly encountered. Phospholipids routinely impact assays as it can be very difficult to separate them chromatographically from drug molecules. They tend not to elute as discrete peaks during the LC analysis, therefore co-eluting with analytes throughout the chromatogram. Sample preparation can be used to reduce these matrix effects through selective removal of phospholipids from the sample, prior to analysis. The efficiency of phospholipid removal can vary significantly depending on the sample preparation technique used.

This presentation aims to compare the effectiveness of different sample preparation techniques (protein precipitation, supported liquid extraction, and SPE (non-polar and mixed-mode)) for removal of phospholipids from human plasma samples.

## Experimental Procedure

Human plasma samples prepared using each of the sample preparation techniques listed were analyzed by LC-MS. Selected phospholipid ions were monitored in order to investigate the amount of phospholipid remaining in the sample after clean-up, and hence evaluate the effectiveness of each technique in removal of phospholipids from plasma samples.

Post-column infusion using propranolol was used to identify the areas of the chromatogram most affected by matrix effects due to the phospholipids.

### Sample Preparation

100 µL of pooled human plasma was extracted using standard methodology for each technique (details below). Following sample preparation, the extracts were evaporated to dryness, and reconstituted in appropriate mobile phase for subsequent analysis. All sample preparation procedures were carried out using 96-well format sample preparation plates.

#### a) Protein Precipitation

ISOLUTE® Array PPT+ (p/n 120-2040-RP)

1. Add 300 µL acetonitrile to each well.
2. Add 100 µL plasma to each well.
3. Allow to stand for 2 minutes.
4. Apply vacuum at -20 "Hg and collect filtrate.

#### b) Supported Liquid Extraction (SLE)

ISOLUTE SLE+ 200mg (p/n 820-0200-P01)

1. Mix 100 µL of plasma with 100 µL H<sub>2</sub>O and apply to each well.
2. Apply a short pulse of vacuum to initiate flow.

3. Allow to stand for 5 minutes, and dry with a short pulse of vacuum.
4. Apply MTBE (1 mL) and elute under gravity for 5 minutes. Apply vacuum (2 minutes) to completely remove extraction solvent.

**c) Non-polar SPE (silica based sorbent)**

ISOLUTE Array C2, C8, C18 25 mg/1 mL (320-0025-RP, 290-0025-RP 220-0025-RP respectively)

1. Condition each well with methanol (1 mL).
2. Equilibrate with ammonium acetate buffer (0.05 M, pH 6, 250 µL).
3. Load plasma sample (100 µL diluted 1:1 (v/v) with ammonium acetate buffer (0.05 M, pH 6)).
4. Wash with buffer:methanol (95:5, v/v, 250 µL).
5. Elute with methanol / 1M ammonium acetate 99.5:0.5 v:v (250 µL).

**d) Non-polar SPE (polymer-based sorbent)**

EVOLUTE™ ABN Array 25 mg/1 mL (p/n 600-0025-RP)

1. Condition each well with methanol (1 mL).
2. Equilibrate with 0.1 % formic acid (1 mL).
3. Load plasma sample (100 µL diluted 1:3 (v/v) with 1% formic acid).
4. Wash with water:methanol (95:5, v/v, 1 mL).
5. Elute with methanol (500 µL).

**e) Mixed-mode SPE (silica-based sorbent)**

ISOLUTE Array HCX 25 mg/1 mL (p/n 902-0025-RP)

1. Condition each well with methanol (1 mL) .
2. Equilibrate with ammonium acetate buffer (0.05 M, pH 6, 250 µL).
3. Load plasma sample: (100 µL diluted 1:1 (v/v) with ammonium acetate buffer (0.05 M, pH 6)).
4. Wash with:
  - (i) ammonium acetate buffer (0.05 M, pH 6, 250 µL),
  - (ii) acetic acid (1 M, 250 µL), dry sorbent with 30sec pulse of vacuum,
  - (iii) methanol (250 µL).
5. Elute with methanol:NH<sub>4</sub>OH (95:5, v/v, 2 x 100 µL).

**HPLC Conditions**

**Instrument:** Waters 2795 Separations module.

**Column:** Zorbax Eclipse XDB-C18. (50 x 2.1mm, 3.5 µL), Agilent

**Guard Column:** Zorbax Eclipse XDB-C8. (12.5 x 2.1, 5 µL)

**Injection volume:** 10 µL. The entire column effluent was directed into the MS.

**Table 1.** Gradient

Time	A = 0.1% Formic in water	B = acetonitrile	Flow rate mL/min
0	70	30	0.25
6	10	90	0.25
8	10	90	0.25
8.1	70	30	0.25

**MS Conditions**

Instrument: Waters Quattro Ultima Pt triple quadrupole MS equipped with an electrospray source.

Source Temp: 100°C

Desolvation Temp: 350°C

Post column infusion: analyte (propranolol) concentration of 1 ng/µL, infused at a rate of 5 µL/min.

**Table 2. SIR Conditions**

Analyte	SIR Mass	Dwell Time (s)	Cone Voltage (V)
Propranolol (PCI Analyte)	260.1	0.1	55
Phospholipid ions	496	0.1	55
	520	0.1	55
	522	0.1	55
	524	0.1	55
	760	0.1	55
	786	0.1	55
	806	0.1	55

## Results

### (i) Matrix effects due to phospholipids

The phospholipids investigated in this presentation were found to have a significant suppression effect on propranolol. This is most apparent in the 5.5–8.0 minute range, shown on the mass chromatograms below (**Figure 1**). The phospholipid peaks in this area are shown on the TIC trace (top, consisting of the ions listed in **Table 2**). The mass chromatogram for propranolol (bottom,  $m/z$  260.1) shows significant suppression correlating to the retention times observed for the phospholipid peaks.




Figure 1. TIC for selected phospholipid ions, with corresponding trace for propranolol (post column infusion,  $m/z$  260.1). Sample: human plasma prepared using protein precipitation (ISOLUTE PPT+) as described above.

## (ii) Comparison of effectiveness of different sample preparation techniques in removal of phospholipids from human plasma samples

**Figure 2** shows the TIC for residual phospholipids in human plasma samples prepared using each of the sample preparation techniques listed. Significant variation in the effectiveness of the different techniques for removal of phospholipids can be observed. This is shown graphically in **Figure 3**.




Figure 2. TIC of residual selected phospholipid ions in human plasma prepared using the various sample preparation techniques




Figure 3. Area count (TIC, 5.5-8.0 minutes) for residual selected phospholipid ions in human plasma prepared using the various sample preparation techniques

## Discussion

1. ISOLUTE PPT+ was least effective in removal of phospholipids, as this technique is based on protein removal only.
2. ISOLUTE SLE+ gave the cleanest extract (most effective phospholipid removal). This probably due to the low solubility of phospholipids in the extraction solvent (MTBE) used.
3. Effectiveness of phospholipid removal using non-polar SPE was inversely proportional to the hydrophobic character of the sorbent. Thus EVOLUTE ABN extracts (polymer based sorbent with high hydrophobic character) contained relatively high amounts of phospholipids, whereas ISOLUTE C2, with very low hydrophobic character, gave very clean extracts.
4. ISOLUTE HCX, a mixed-mode sorbent, gave relatively high phospholipid content, despite the rigorous interference elution regime possible. This is believed to be due to interactions of the zwitterionic phospholipids with both the non-polar and cation exchange functional groups. Further work will investigate this result.

## Overall Conclusions

Sample preparation techniques can be ranked as follows for their effectiveness in removal of phospholipids from human plasma samples (from most effective to least effective).

- Supported liquid extraction
  - ISOLUTE SLE+
- Non-polar SPE (silica)
  - ISOLUTE C2
  - ISOLUTE C8
  - ISOLUTE C18
- Mixed-mode SPE
  - ISOLUTE HCX
- Non-polar SPE (polymer)
  - EVOLUTE ABN
- Protein precipitation
  - ISOLUTE PPT+

© 2006 Argonaut Technologies, now a Biotage company. All rights reserved. ISOLUTE is a registered trademark and EVOLUTE is a trademark of Argonaut Technologies, now a Biotage company.

**[www.biotage.com](http://www.biotage.com)**

### United States and Canada

Tel: +1 434 979 2319  
Toll-Free: +1 800 446 4752  
[ordermailbox@biotage.com](mailto:ordermailbox@biotage.com)

### United Kingdom, EIRE

Biotage  
Tel: +44 1992 501535  
[order@eu.biotage.com](mailto:order@eu.biotage.com)

### Sweden

Biotage  
Tel: +46 18 56 59 00  
[order@eu.biotage.com](mailto:order@eu.biotage.com)

### Japan

Biotage  
Tel: +81 422 281233  
[order@biotage.co.jp](mailto:order@biotage.co.jp)