

# **Automating TLC to Flash Purification Gradient Methods Containing Polar Solvents\***

**Sjaan Armentrout, Jack Liu and Peter Rahn**

**Biotage**

**Discovery Chemistry Group**

**1725 Discovery Drive**

**Charlottesville, VA 22911**


**\* Patent Pending**



# Abstract

- **Biotage SP1 and SP4 systems automatically produce appropriate optimized Flash gradient profiles from TLC data**
- **These automated gradient methods are applicable for a wide range of solvent systems including those containing acetonitrile or alcohols**
- **Examples of reductive amination products utilizing TLC data to produce successful Flash gradient purification separations are presented**

# TLC Data Gathering



- **New Biotage Flash systems automatically utilize TLC data to generate optimized Flash purification conditions**
- **Eliminates need to optimize TLC separation**
- **Chemist empirically determines suitable solvent composition to separate target compound from nearest impurity**

**Calculate  $R_f = b/a$**

**Identify Solvent and ratio X:Y**


# New Automated Biotage Flash Systems



- **Chemist Inputs the following TLC Data**
  - **R<sub>f</sub> of selected compound**
  - **R<sub>f</sub> of additional compounds**
    - If R<sub>f</sub> of more than one compound is entered system determine sample load or cartridge size)
  - **Solvent identity and ratio**
- **Automatic Flash System Setup\***
  - **Identifies weaker solvent**
  - **Produces appropriate gradient profile including initial hold, gradient and final hold**
  - **System suggests altering solvent composition from TLC concentration if selected solvents exceed a specified solvent strength threshold**

\* Patent Pending

# Automatic Optimization of Purification Conditions for Any Component In a Mixture




**TLC Conditions**

**DCM:MeOH  
98:2**


- **Algorithm works for Non-Polar and Polar Solvents**
  - MeOH was required to separate these polar compounds
  - If TLC solvents include polar solvents such as acetonitrile, methanol, ethanol, isopropanol etc. the system automatically notifies chemist to adjust the solvent mixture to an appropriate value
- **Rf factors between 0.05 and 0.9 are suitable**
  - Eliminates the need to optimize Rf between 0.15 and 0.3
- **Automated Sp1 System optimizes gradient conditions to purify identified product**
  - Includes initial, gradient and final hold conditions
  - Suggests adjustments in mobile phase composition
  - Suggests column load or column size if more than one Rf is entered

# Optimized Flash Gradient Profile to Purify 4-Iodoaniline




- **R<sub>f</sub> of 4-Iodoaniline= 0.54** entered into SP1 system and **DCM:MeOH** solvent system identified and **Flash 12+ cartridge** specified
  - **Automatic gradient profile generated:**
    - Initial conditions = 18 mL at 0% B
    - Linear gradient = 120 mL 0-46%B
    - Final conditions = 36 mL at 46% B
- Resolves 4-Iodoaniline and 2-Benzyl-2-DMBP**
- Solvent strength wasn't sufficient to elute 4-DMAAP- required only 1 automatic extension ( - - - - - )**

# Optimized Flash Gradient Profile for 2-Benzyl-2-DMBP



- **R<sub>f</sub>=0.40** for 2-Benzyl-2-DMBP entered into SP1 system
- **Resulting gradient:**
  - Initial conditions = 18 mL at 0% B
  - Linear gradient = 120 mL 0-60%B
  - Final conditions = 36 mL at 60% B
- **Achieved excellent resolution of first 4-Iodoaniline and 2-DMBP**
- **Solvent strength still not sufficient to elute 4-DMAAP and an automatic gradient extend was initiated (-----)**


# Optimized Flash Gradient Profile to Purify 4-DMAAP



- **$R_f=0.11$  for 4-DMAAP (compound 3) entered into SP1 system**
- **Resulting automatic gradient profile:**
  - Initial conditions = 18 mL at 0% B
  - Linear gradient = 120 mL 0-89%B
  - Final conditions = 36 mL at 89% B
- **Good resolution of all three compounds, no gradient extension required**
  - Even with this stronger solvent mixture all the polar compounds remained resolved



# SP1 Flash Purification Results Using Automated Purification Conditions for Polar Compounds



Fr 2/3 Fr 7/8 Fr 10/11 4-1 DMPB DMAAP

**TLC Conditions**

**DCM:MeOH 98:2**

- TLC shows fraction analysis after Automated Flash Purification using DCM:MeOH
- Algorithm works for Rf factors between 0.05 and 0.9
- Automated Sp1 System optimizes gradient conditions to purify identified product
  - Includes initial, gradient and final hold conditions
  - Suggests adjustments in mobile phase composition
  - Suggests column load or column size if more than one Rf is entered
- MeOH was required to separate these polar compounds
  - System automatically notifies user to modify polar solvent ratio
- Examples using this patented technology for the purification of amination products in shown in the next section of this poster

# Synthesis and Purification of Amination Products using Biotage Initiator™ Microwave Synthesis and Biotage SP1™ Flash Purification Systems



- **Biotage Initiator™ 60 Microwave System** was equipped with the following:
  - Compact footprint
  - Touch screen control
  - Four different vial sizes in any order or combination at any time without system modifications for greater flexibility and direct scale-up of milligrams to grams
  - 60-position sample bed
  - Best-in-class safety




- The 2-5 mL reaction vials containing the solid supported reagent, solutions and enclosed magnet to stir reaction mixture during synthesis



- Purifications performed on the compact SP1 system using the new TLC to gradient algorithm
  - Accommodates the Biotage FLASH 12-40 cartridges without external stands
  - Solvent and waste monitoring
  - Variable dual-wavelength detector
  - Auto-continue feature


# Microwave Synthesis and Purification of Reductive Amination Product



**Step 2** Microwave 5 minutes @ 150 °C

**Step 3** Confirm product formed by TLC evaluation

**Step 4** FLASH chromatography with catch and release using Syntage TsOH Samplet and DCM:Methanol gradient



- Starting amine and ketone materials placed in reaction vessel with 20% HOAc in THF as solvent and Syntage Cyanoborohydride reagent

- Reaction mixture heated and stirred in Initiator microwave at 150°C for 5 minutes

- Reaction mixture tested by TLC

- Reaction mixture solution containing THF directly added to 12+M cartridge containing a TsOH Samplet

- From TLC input Rf data for compound of interest

- Product automatically and successfully purified on SP1 Flash system using catch and release technique

# Syntage TsOH Samplet Catch and Release Procedure

Biotage Si 12+M cartridge containing TsOH  
Catch and Release Scavenger Samplet

## Step 1

Add reductive amination reaction mixture in THF to Samplet/Cartridge



## Step 2

Capture product and wash out THF with DCM

Discard THF washed from column with DCM

Reaction product (amine) bound to TsOH Samplet


## Step 3

Release reaction product and impurities from Syntage Samplet with 0.5 mL 2.0M  $\text{NH}_3$  in MeOH

## Step 4

Purify selected compounds with specified DCM:MeOH gradient


# TLC and Purification Data for Reductive Amination Product with DCM and MeOH



## HPLC Conditions


Column: C<sub>18</sub> 5μm 4.6 x 50 mm  
Mobile Phase: A. H<sub>2</sub>O:ACN (95:5) + 0.1% TFA  
B. ACN:H<sub>2</sub>O (98:2) + 0.1% TFA  
Gradient: 0-100% B in 8 min, hold 2 min  
Flow Rate: 1.5 mL/min  
Detection: 254nm

## TLC analysis of crude reaction mixture



- **TLC solvent system was DCM:MeOH 97:3**
  - Product Rf 0.34
- **Gradient profile automatically generated by Flash system**
  - Solvent A: DCM
  - Solvent B: DCM:MeOH (94:6)
  - Initial hold for 18 mL at 0% B
  - 120 mL 0-66% B linear gradient
  - Final hold for 36 mL at 66% B

# Reductive Amination Product Purification Results with DCM and MeOH



- **Purification with Biotage Flash SP1 system using Automated DCM/MeOH gradient**


- Flash conditions: Flash 12+M cartridge loaded with 110 mg crude

Solvent A: DCM

Solvent B: DCM:MeOH (94:6)

Gradient: Initial hold for 18 mL at 0% B  
120 mL 0-66% B linear gradient  
Final hold for 36 mL at 66% B

- **Total Purification Time 15 minutes**



**HPLC Conditions**

Column: C8 4.6 x 50 mm

Mobile Phase A: H<sub>2</sub>O:ACN (95:5) + 0.1% TFA

Mobile Phase B: ACN:H<sub>2</sub>O (98:2) + 0.1% TFA

Gradient: 5-100% B in 8 minutes, hold 2 min

Flow Rate: 1.5 mL/min

Detection: 254nm

# Advantages of New Flash Automated Gradient Technique

- **Accelerates Solvent System Selection**
  - Eliminates guessing the solvent mixture that produces both an acceptable separation and a  $R_f$  within a narrow band
  - Eliminates need to perform multiple TLC separations
  - Technique applicable when  $R_f$  ranges from 0.05 to 0.9
  - New Algorithms applicable for all solvents
- **Provides a quick, easy solution for alcohol containing gradients**
  - Good for very polar compounds that require alcohol
  - System automatically adjusts conditions using algorithms when solvents exceed a predetermined solvent strength value with no user intervention
- **Eliminates remixing solvents to elute different compounds**
  - User optimizes separation of mixture components by inputting different  $R_f$  values
- **Accelerates Clean-up and Purification**
  - Biotage catch and release scavenger Samplets eliminate time-consuming clean-up of reaction mixture prior to purification
  - Catch and release technique ensures polar solvents do not cause co-elution or early breakthrough of desired compounds