


Microwave-Assisted Sulfamide Synthesis

Shahnaz Ghassemi^a, Kristin Fuchs

Biotage, Discovery Chemistry Group 1725 Discovery Drive Charlottesville, VA 22911

Introduction

In the area of combinatorial library synthesis for medicinal agents there is a constant need for new methodologies [1]. The sulfamide compounds are noted for their broad and potent antibacterial activity [2-3]. The unsymmetric sulfamides appear to be more potent as protease inhibitors than the symmetric analogues due to the flipped conformation that occurs during binding [4]. Unfortunately, most syntheses focus on symmetric sulfamides. The few methods available for unsymmetrical compounds rely on low-yielding synthetic steps that are neither general nor selective [1,5]. A novel transition-metal-catalyzed process for making unsymmetric sulfamides that was recently reported has several limitations, especially with ortho-isomers [1]. Even though other available methods report high yields, they either require reagents that are not readily accessible or they focus on specific structures rather than a general procedure [6]. Winum and co-workers reported a novel sulfamoylating reagent used in the synthesis of sulfamides [8]. However, our study showed that using the sulfamoylating reagent added additional steps and resulted in slow, low-yielding reactions. In an effort to find a fast and general method for sulfamide synthesis we found that microwave heating facilitates the synthesis of sulfamides. This was accomplished in one-pot reaction by a stepwise addition of CSI to tertbutanol at 0 °C to form the N-(tert-butoxycarbonyl) sulfamoyl chloride intermediate 2 (scheme 1). Anilines or amines were added the reaction mixture was heated using microwave heating at 80 °C for five minutes. The resulting products were isolated using normal-phase flash chromatography with a good yield (table 1). Microwave synthesis provided great improvements in increasing product yield and decreasing reaction time [7].



The microwave assisted Mitsunobu reaction was used for alkylation of Boc-sulfamides with different alcohols (Scheme 2) [9]. The reaction time depended on the structure of alcohols. For example, microwave irradiation of benzyl alcohol mixture with Boc-sulfamides, triphenyl-phosphine and diethyl azodicarboxylate (DEAD) in THF provided N-alkylated products in four minutes at 80 °C. In the case of 2-pyridinepropanol (compounds 10, 12, 14), one minute of microwave heating was enough to complete reaction.




Table 1: One-pot microwave-assisted sulfamides synthesis

#	Amines	Products ^a	% Yield ^b	% Purity ^c	(M ⁺ + Na)
3			80	96	309.93
4			55	97	295.94
5			62	93	308.93
6			67	96	334.97
7			86	97	363.93

^aReactions were performed in the Biotage EMRY5™ Liberator microwave system in 2-5 mL reaction vials at 80 °C. ^bYield of isolated product. All products were isolated on the Biotage Sp automated flash chromatography system (Flash 25+ M, 25 x 150, 40-63 μm, 60 Å), using ethyl acetate and hexane gradient. ^cPurity is calculated by HPLC (Waters, C8 4.6 x 50 mm, S-3 120 Å). ^dMass spectroscopy was carried out on a Micromass® ZQ (Waters). ^e¹H NMR data in CDCl₃ was collected on a 500 MHz Bruker spectrometer.

The tert-butoxycarbonyl group removal is generally carried out with trifluoroacetic acid either neat or in combination with CH₂Cl₂ [14]. Since CF₃COOH is volatile, harsh and corrosive, a search for an alternative method of deblocking is ongoing.

Recently was reported that Amberlyst 15, a strong acidic resin, can remove the Boc-protecting group and form salts with the deprotected amines [15]. This method has been used to facilitate the generation and purification of amines. However, this technique requires a long reaction time (12-24 hours). We decided to explore the scope and limitations of deblocking the BOC-group from sulfamides using silica-bonded phenylsulfonic acid, and the effects of microwave heating in altering the reaction time (scheme 3). Boc-sulfamides were treated with Si-TsOH and heated by microwaves at 100 °C. In all the examples, the Boc-protecting group was completely removed within five minutes [16]. Here we report that microwave heating with Si-TsOH significantly shortens the Boc-removal time. The formation of salts between the sulfamide and silica bonded acid depends on substituents on the sulfamide nitrogen (pKa of sulfamides 7-11). The desired products were released from Si-TsOH surface using NH₃/MeOH (scheme 3).




Table 3: Microwave-assisted BOC- deblocking using Si-TsOH

- General procedure for Boc-sulfamides using microwave heating: In a typical experiment, chlorosulfonyl isocyanate (0.24 ml, 2.7 mmol) was added dropwise to a solution of tert-butyl alcohol (0.26 ml, 2.7 mmol) in anhydrous dichloromethane (3 ml) in a sealed Pyrex tube under inert gas at 0 °C. Amine (5.5 mmol) was then added and the reaction was heated in a microwave cavity for 5 minutes at 80 °C. The reaction mixture was added to a Samplet™ cartridge and purified by flash chromatography.
- General procedure for microwave assisted Boc-sulfamides cleavage with Si-TsOH: Method A: Silica-bound p-toluenesulfonic acid (1.26g, 0.96 mmol) was added to the Boc-protected sulfamide (0.32 mmol) in 1:1 acetonitrile: DCM (4ml). The reaction was heated to 100 °C in a microwave cavity for 5 minutes. The reaction mixture was then loaded onto a silica column. Using the following conditions on flash chromatography yielded the desired compound.

#	Products ^a	% Yield ^b	M ⁺ + Na ^c
15		90	209.11
16		75	195.03
17		75	209.05
18		61	235.63
19		91	264.02
20		71	325.38
21		92	353.99
22		89	354.42

^aReactions were performed in the Biotage EMRY5™ Liberator microwave system in 2-5 mL reaction vials. ^bYield of isolated product. ^cMass spectroscopy was carried out on a Micromass® ZQ (Waters). ^d¹H NMR data was collected on a 500 MHz Bruker spectrometer.

Conclusion

A general microwave assisted reaction in preparation of unsymmetric Boc-sulfamides is demonstrated. Also, an alternative method of Boc-removal from sulfamide was introduced using Si-TsOH in conjunction with microwave heating. Boc-deblocked sulfamides were captured by Si-TsOH, depending on their pKa. The captured sulfamides were released from surface of Si-TsOH by using NH₃ in MeOH, followed by quick flash purification. This new method of microwave-assisted, Boc-cleavage group from sulfamides facilitates the preparation and purification of unsymmetric sulfamides.

References:

- Alcaraz, L.; Bemmion, C.; Morris, J.; Meghani, P.; Thom, S. *Org. Lett.* (2004), 0, A-D
- Nishino, Y.; Komarasaki, T.; Yuasa, T.; Kakimura, M.; Izuno, K.; Kobayashi, T.; Masui, Y.; Fujie, S.; Gotoh, T.; Hajima, M.; Takahira, M.; Kataoka, T.; Okuyama, A. *Organic Process Research and Development*, (2003) 7, 846-850.
- Markgren, P.; Schaal, W.; Hamalainen, M.; Karlen, A.; Hallberg, A.; Samuelsson, B.; Danielson, H. J. *Med. Chem.* (2002), 45, 5430-5439.
- Schaal, W.; Karlsson, A.; Ahlson, G.; Lindberg, J.; Andersson, H.; Danielson, H.; Classon, B.; Uнге, T.; Samuelsson, B.; Hulén, J.; Hallberg, A.; Karlen, A. *J. Med. Chem.* (2001) 44, 153-169
- Hof, F.; Irvine, P.; Johnson, D.; Robek, J. *Org. Lett.* (2001), 3, 4247-4249
- Masui, T.; Kabaki, M.; Watanabe, H.; Kobayashi, T.; Masui, Y. *Organic Process Research and Development*, (2004) 8, 408-410.
- Alcaraz, L.; Bemmion, C.; Morris, J.; Meghani, P.; Thom, S. *Org. Lett.* (2004), 0, A-D



EMRY5™ Liberator Microwave System



SP1 Automated Purification System

Biotage