

Method Development Strategies Using Polymer-based SPE for the Analysis of Peptides Prior to UPLC-MS/MS Analysis

Helen Lodder¹, Lee Williams¹, Rhys Jones¹, Alan Edgington¹, Adam Senior¹, Geoff Davies¹, Steve Jordan¹, Claire Desbrow¹, Paul Roberts¹, Stephanie Marin², Dan Menasco² & Elena Gairloch²

¹Biotage GB Limited, Distribution Way, Dyffryn Business Park, Ystrad Mynach, Cardiff, CF82 7TS, UK

²Biotage, 10430 Harris Oaks Blvd., Suite C, Charlotte North Carolina 28269, USA

Introduction

Two peptides, produced in the hypothalamus, oxytocin a 9-residue neuropeptide associated with parturition (childbirth) and lactation and vasopressin a peptide hormone associated with fluid retention and vasoconstriction were analyzed, **Figure 1**. This poster investigates the selectivity of polymer-based solid phase extraction sorbents and rapid extraction of clinical matrices, matrix interference reduction and improved analyte detection.




Figure 1. Structure of vasopressin (left) and oxytocin (right).

Experimental

Reagents

Oxytocin and vasopressin were obtained from Phoenix Peptide. Ammonium hydroxide, ammonium acetate, formic acid and LC/MS grade solvents were purchased from Sigma-Aldrich Chemical Co. (Poole, UK). Water (18.2 MΩ.cm) was drawn fresh daily from a Direct-Q 5 water purifier (Merck Millipore, Watford, UK). Human serum was purchased from Golden West Biologicals Inc. (Ca, USA.). All materials for gel electrophoresis were obtained from ThermoFisher Scientific.

Sample Preparation

Extractions were performed using the 30 mg 96-well plate format.

SPE Optimization: SPE investigated polymer-based non-polar retention with corresponding mixed-mode strong and weak cation exchange mechanisms using the EVOLUTE[®] EXPRESS sorbent family. The final SPE protocols are detailed in **Table 1**. All solvent/buffer steps were optimized to 1 mL while elution solvent volumes were minimized to 200 µL.

Oxytocin and vasopressin: 200 µL of serum spiked at 50 ng/mL.

Table 1. Optimized SPE procedures.

Step	EVOLUTE [®] EXPRESS [®] ABN	EVOLUTE [®] EXPRESS [®] CX	EVOLUTE [®] EXPRESS [®] WCX
Condition	MeOH	MeOH	MeOH
Equilibration	0.1% Formic acid	50 mM NH ₄ OAc pH6	H ₂ O
Pre-treatment	1:1 1% Formic acid	1:1 0.1% NH ₄ OH	1:1 1% NH ₄ OH
Sample load	400 µL	400 µL	400 µL
Wash 1	0.1% Formic acid	0.1% NH ₄ OH	1% NH ₄ OH
Wash 2	-	90:10 H ₂ O: MeOH	90:10 H ₂ O: MeOH
Elution	5% FA 80/20 H ₂ O/MeCN	2% NH ₄ OH in MeOH	2% FA 30:70 H ₂ O/MeCN

Post extraction: Extracts were either evaporated to dryness using an SPE Dry unit at 40 °C and reconstituted in 200 µL 0.1% FA 90/10 aq/MeCN or direct injected.

Gel Electrophoresis Procedure

Details available on request.

UPLC Conditions

Instrument: Waters ACQUITY I-Class UPLC equipped with a 15 µL flow-through needle (Waters Assoc., Milford, MA, USA)

Column: ACE C18-300: 50 mm x 2.1 mm id, 1.7 µm, (ACT, UK)

Mobile Phase: **A:** 0.1% FA aq; **B:** 0.1% FA MeCN

Flow Rate: 0.4 mL/min

Gradient: Linear ramp 10-33% B over 1.5 min; resume initial starting conditions at 1.6 min

Column Temperature: 40 °C

Injection Volume: 10 µL

Mass Spectrometry

Instrument: Xevo TQ-S triple quadrupole mass spectrometer (Waters Assoc., Manchester, UK) equipped with an electrospray interface.

Positive ions were acquired in the multiple reaction monitoring (MRM) mode as shown in **Table 2**.

Desolvation Temperature: 500 °C; **Ion Source Temperature:** 150 °C

Collision Gas Pressure: 3.6 x 10⁻³ mbar

Table 2. MS conditions.

Compound	MRM Transition	Cone Voltage (V)	Collision Energy (eV)
Oxytocin (Qual)	1007.3 > 285.2	30	25
Oxytocin (Quant)	504.2 > 285.2	30	15
Vasopressin (Qual)	542.8 > 157.2	30	10
Vasopressin (Quant)	542.8 > 328.2	30	15

Results

LC-MS/MS Optimization

Oxytocin and vasopressin demonstrated precursor ions at 504.2 and 542.8 m/z, respectively, corresponding to [M + 2H]²⁺ multiple charging. Chromatographic performance was compared using various phase chemistries as shown in **Figure 2**. Optimum chromatography was obtained using the 50 mm ACE C18-300 Å phase.




Figure 2. Column chromatography comparison.

Extraction Optimization

SPE strategies focused on extraction pH, wash and elution solvent combinations for optimal extraction recoveries and cleanliness. The CX chemistry demonstrated minimal extraction of vasopressin and around 60% recoveries of oxytocin. Further method development demonstrated no improvement. ABN and WCX chemistries delivered greater than 60% recoveries for both analytes as shown in **Figure 3**.




Figure 3. Extraction recovery profiles for SPE chemistry comparison.

Pre-treatment and wash solvent optimization for ABN and WCX chemistries are demonstrated in **Figure 4**.




Figure 4. Extraction recovery profile for pre-treatment and wash solvent optimization using ABN (left) and WCX (right).

Final modification to improve extract cleanliness involved pH modification and elution in high aq proportions. **Figure 5** demonstrates acid comparison and maximum aq content for each ABN and WCX chemistries.




Figure 5. Extraction recovery profile for elution optimization using ABN (left) and WCX (right).

The presence of water wettable components allows the possibility to eliminate plate conditioning steps prior to sample loading resulting in a simple load-wash-elute procedure. **Figure 6** demonstrates final extraction recoveries for full method (documented in **Table 1**) and L-W-E procedures combining traditional evaporation/reconstitution and direct injection using EVOLUTE[®] EXPRESS ABN.




Figure 6. Extraction recovery profile for optimized extraction conditions using EVOLUTE[®] EXPRESS ABN.

Final extraction protocols were investigated for protein removal and phospholipid content. **Figure 7** demonstrate gel electrophoresis profile of ABN extracts using various elution volumes.

Figure 7. ABN Gel electrophoresis protein profile.




Figure 8. demonstrates the total ion chromatograms (TICs) of typical phospholipid MRMs obtained comparing ABN and WCX extracts.





Figure 8. Phospholipid MRM TICs for final serum extracts for each chemistry.

Final extract cleanliness was further investigated using post-column infusion (PCI) experiments. Oxytocin and vasopressin were infused post-column, directly prior to the ESI source to allow easy identification of regions of suppression throughout the assay, as shown in **Figure 9**.

Figure 9. PCI baselines of the optimized ABN extraction protocol: blank solvent (red), extracted serum using full SPE (purple) and direct inject procedure (green).



EVOLUTE[®] EXPRESS ABN demonstrated overall better performance. Spiked calibration curves constructed from 25-1000 pg/mL returned excellent coefficients of determination (r²) and sensitivity. Typical calibration lines are demonstrated in **Figure 10**.




Figure 10. Serum calibration curves extracted using EVOLUTE[®] EXPRESS ABN.

Conclusion

- EVOLUTE[®] EXPRESS WCX provided good recoveries. However, better cleanliness was obtained using the optimized elution combinations with EVOLUTE[®] EXPRESS ABN.
- SPE optimization resulted in recoveries greater than 70% while demonstrating good removal of matrix components in the form of proteins and phospholipids.
- It was also possible to eliminate evaporation and move to direct inject using the selected chromatographic mobile phases, however, there was an increase in regions of suppression observed.
- Optimized SPE protocol still showed excellent protein removal whilst maintaining good peptide recoveries.