

# The Application of MAOS to Medicinal Chemistry and Drug Discovery

David Lonergan, PhD  
Director, Medicinal Chemistry  
Nanosyn, Inc.

# Why Use MAOS?

- Accelerated reaction kinetics
- Increased conversion of starting materials to products
- Reaction mixtures are often cleaner than conventional synthesis
- Small to large scalability
- Auto-sampler allows for overnight runs
- Extremely easy to operate

# Why is MOAS Important for Nanosyn?

Nanosyn is a contract research organization (CRO)

- Multiple projects with multiple customers at any given time
- Chemistry constantly changes from one day to the next
- Non-trivial synthesis: heterocyclic ring formation, C-C and C-N bond formation, multi-step (10-12 steps)
- Synthetic method development needs to be completed in a very short time frame—typically 1-2 days

## How Often is MAOS Applied at Nanosyn?

- MAOS is used almost daily by every chemist at Nanosyn
- During the last 4 years, chemists used MAOS an average of 1.3 times per day
- The microwave is typically the first place we start when investigating a new chemical synthesis or method development

# What Types of Chemistry Have Been Applied to MAOS at Nanosyn?

- Palladium catalyzed reactions  
*Suzuki, Stille, Buchwald, Sonogashira*
- Nucleophilic aromatic substitution reactions
- Reductive amination
- Amide bond formation
- Heterocycle formation

We have applied MAOS to almost every type of reaction possible

# What Types of Chemistry Have Been Applied to MAOS at Nanosyn?

- Multi-step scaffold synthesis
- Multi-step single compound synthesis
- Mg to multi-gram synthesis
- Predominantly C-C and C-hetero-atom and heterocyclic chemistry


*Challenging, diverse chemistry*

# Palladium Catalyzed Reactions

*Suzuki* type reactions:

- One of most important class of contemporary carbon-carbon bond forming reactions
- Require boronic acid and aryl halide or aryl triflate coupling partner
- Not all boronic acids are commercially available
- Conventional conditions require heating at  $\sim 100^{\circ}\text{C}$  for 24 hours


# Palladium Catalyzed Reactions



# How does MAOS Help In Applying the *Suzuki* to MedChem?

The boronic acid is not commercially available


- Conventional boronic acid synthesis:



Drawbacks:

- Ester functionality makes  $n\text{-BuLi}$  use difficult

# MOAS Boronic Ester Synthesis




- Fast and convenient boronic ester synthesis
- Tolerates many functionalities
- Works well with  $Br$ -,  $I$ -, and  $OTf$  derivatives

# MAOS to Provide Access to Unaccessible Substrates


- 2-Substituted pyridyls:

2-pyridyl boronic acids are not commercial available due to proto-deboronation



# Simple Method for Preparing 2-Pyrdyl Substrates

- Stannane in the 2-position is stable




## *Buchwald/Hartwig* Reaction

- Works very well for simple substrates
- Complicated substrates require tedious and difficult optimization
- Very sensitive to temperature, catalyst, base solvent and ligand

**MAOS can be very useful for *Buchwald optimization***

- Multiple conditions can be easily performed
- Excellent control of temperature

# Buchwald/Hartwig Reaction



- Works well with simple substrates
- Heteroaryl amines and halides pose issues
- MAOS can work well when conventional Buchwald fails

# C-N Forming Reactions


## *Reductive amination:*

- Excellent method for alkylating amines
- Major drawback for multiple/parallel reactions is work up

**Combination of MAOS and silica bound cyanoborohydride**

# Reductive Amination


- Fast
- Easy to work-up; centrifuge, decant, evaporate



# MAOS Applied to *Simple* Chemistry


## Amide bond formation

- *Supposed to be simple*
- What about cases with extremely deactivated anilines?



# Amide Bond Formation

- Conventional heating will give only traces of product



# Conclusions

- MAOS has wide applicability to industry oriented medicinal chemistry
- MAOS has changed how Nanosyn's approach synthetic and medicinal chemistry
- MAOS has become indispensable to our method development
- Combining MAOS with excellent synthetic and medicinal chemists is a true formula for success in drug discovery